










Height measuring instruments

	 seca 242 CE 0123	 seca 240 CE 0123	 seca Leicester CE	 seca 217 new CE 0123	 seca 225 CE 0123	 seca 222 CE 0123
	Electronic measuring rod for wall mounting, wireless data transmission onto separate display.	Mechanical precision measuring rod for the wall.	Stand-alone stadiometer for measuring children and teenagers, easy to set up and handy to transport.	Stadiometer for mobile height measuring.	Mobile measuring unit on rollers, for convenient transport	Telescopic measuring rod for the wall with extra-wide measuring slides and heel stop.
Technical Data						
Measuring range	62 – 210 cm / 25 – 82.6"	60 – 210 cm / 25 – 82.6"	20 – 207 cm / 8 – 81.5"	20 – 205 cm / 8 – 81"	6 – 230 cm / 2 – 90"	6 – 230 cm / 2 – 90"
Graduation	1 mm / 1/8"	1 mm / 1/8"	1 mm / 1/8"	1 mm / 1/8"	1 mm / 1/8"	1 mm / 1/8"
Dimensions (WxHxD)	300 x 2,200 x 250 mm	300 x 2,200 x 250 mm	356 x 2,116 x 416 mm	328 x 2,145 x 574 mm	310 x 1,400 x 410 mm	120 x 1,322 x 34 mm
Dimensions/display (WxHxD)	198 x 55 x 157 mm					
Weight	2.5 kg	3 kg	3 kg	3.6 kg	7 kg	1.3 kg
Power supply	Batteries					
Unit price	£ 429	£ 179	£ 55	£ 145	£ 239	£ 125

Multifunction and wheelchair scales

	 seca 645 CE 0109	 seca 685 CE 0109	 seca 677 CE 0109	 seca 665 CE 0109	 seca 657 CE 0109	 seca 675 CE 0109
	Large handrail scales for multifunctional use, high load-bearing capacity and BMI function.	Electronic multifunction scales with stable handrail and high load-bearing capacity, mobile with integrated castors.	Electronic wheelchair scales with large platform, robust handrail and transport castors.	Electronic wheelchair scales with extremely flat step-on surface and ramp, collapsible and mobile.	Electronic scales with large platform for trolleys or stretchers on wheels, with integrated transport castors.	Electronic platform scales with high load-bearing capacity and castors for transport.
Technical Data						
Capacity	300 kg	300 kg	300 kg	300 kg	300 kg	300 kg
Graduation	100 g < 200 kg > 200 g	100 g < 200 kg > 200 g	100 g < 200 kg > 200 g	100 g < 200 kg > 200 g	100 < 200 kg > 200 g	100 g < 200 kg > 200 g
Dimensions (WxHxD)	750 x 1,140 x 750 mm	910 x 1,100 x 965 mm	951 x 1,107 x 965 mm	1,102 x 911 x 1,150 mm	893 x 87 x 1,470 mm	890 x 65 x 965 mm
Dimensions/platform (WxHxD)	600 x 45 x 600 mm	800 x 46 x 850 mm	800 x 51 x 965 mm	800 x 46 x 965 mm	800 x 51 x 1,470 mm	800 x 46 x 965 mm
Weight	24 kg	59 kg	34.5 kg	32.3 kg	43.6 kg	25.5 kg
Power supply	Mains adapter / Batteries	Mains adapter	Mains adapter	Mains adapter / Batteries	Mains adapter	Mains adapter
CE approval class						
Unit price	£ 1,195	£ 1,895	£ 1,495	£ 1,299	£ 1,595	£ 1,095

				
seca 202	seca 216	seca 206	seca 203 new	seca 201 new
Telescopic measuring rod for the wall.	Mechanical stadiometer. Headpiece with locking screw for accuracy.	Roll-up measuring tape with wall attachment.	Ergonomic, smoothly extractable circumference measuring tape with WHR calculator to determine body fat distribution.	Ergonomic, smoothly extractable measuring tape for determining body circumference.
6 – 200 cm / 2 – 78"	150 cm / 60" selectable from 10 – 230 cm	0 – 220 cm / 0 – 87"	15 – 205 cm / 6 – 80"	15 – 205 cm / 6 – 80"
1 mm / 1/8"	1 mm / 1/8"	1 mm or 1/8"	1 mm or 1/8"	1 mm or 1/8"
55 x 1,322 x 34 mm	120 x 1,500 x 217 mm	125 x 125 x 170 mm		
900 g	2 kg	140 g	about 75 g	about 55 g
£ 69	£ 59	£ 17	£ 60 (Pk of 10)	£ 30 (Pk of 10)

Accessories for multifunction and wheelchair scales

			
seca 472	seca 223	seca 223	seca 470
Stand for display elements of seca scales with cable remote displays.	Telescopic measuring rod for the seca 645 handrail scales.		Drive-up ramp as accessory for seca 675.
Technical Data			
Measuring range		15 – 230 cm / 6 – 90"	
Graduation		1 mm / 1/8"	
Dimensions (WxHxD)	420 x 1,030 x 420 mm	120 x 1,322 x 34 mm	805 x 29 x 250 mm
Weight	7.2 kg	1.3 kg	1.6 kg
Unit price	£ 79	£ 129	£ 75