

LABCOLD MODEL RLDF18041

505 Litre Pharmacy Refrigerator



- Conforms to RPSGB guidelines
- Conforms in full with Department of Health "Green Book" requirements.
- Digital temperature display
- Integral high/low recording
- Fan assisted air circulation
- Door lock
- CFC, HCFC and ammonia free
- Automatic defrost
- Accurate temperature control +2°C to +8°C
- Audio alarm with mute facility

Large vaccine/pharmacy refrigerator designed for the volume storage of medicines and vaccines. This model, exclusively developed by Labcold is complete with all the technology required to ensure expensive stock is kept in optimum condition including features such as integral temperature display, alarm and forced air cooling.

Completely compliant with all the regulations specified by both the Department of Health and the RPSGB this refrigerator is the professional choice for safe temperature controlled volume vaccine and pharmaceutical storage.

SPECIFICATIONS - RLDF18041	
Configuration	Free standing
Capacity	505 litres
Exterior size (HxWxD) (mm)	1695 x 770 x 730
Interior Fan	✓
Door Lock	✓
Alarm	✓
Auto Defrost	✓
Shelves	5
Operating temp (nominal)	+2 to +8°C
Lock type	Manual key lock

Supplied By MediSafe Ltd

MediSafe Ltd · Kings Worthy House · London Rd · Kings Worthy · Winchester · SO23 7QA
Tel: 01962 885777 · Fax: 01962 885678 · Email: sales@medisafe.org.uk · Web: www.medisafe.org.uk