

LABCOLD MODEL RPDF0012D

Portable Laboratory & Medical Transporter



- Temp range:-5 to +60°C
- Two speed fan with mute facility
- Lightweight
- CFC, HCFC and ammonia free
- Digital set point display

Simple and versatile, the Labcold Digital Portable Laboratory & Medical Transporter can both cool or heat, and is designed for the transporting of samples, vaccines, pharmaceuticals or anything that needs to be moved under temperature controlled conditions.

Complete with an in-car power cable as mains, this transporter can keep contents at the desired temperature whether in transit or in situ.

This product is not suitable for the transport of blood or blood products.

SPECIFICATIONS - RPDF0012D	
Configuration	Electronic Portable
Capacity	12 Litres
Exterior size (HxWxD) (mm)	290 x 415 x 285
Interior Fan	✓
Temperature range	-5 to +60°C
Door Type	Top opening
Compartments	2
Weight	6 Kg

Supplied By MediSafe Ltd